## 3M ESPE

## Ketac™ Silver

- EN** **Silver-Reinforced Glass Ionomer Restorative Material**
- DE** **Silberverstärktes Glasionomer Füllungsmaterial**
- FR** **Matériau d'obturation en verre ionomère à renfort d'argent**
- IT** **Materiale per otturazioni vetro-ionomero rinforzato all'argento**
- ES** **Material de obturación de vidrio ionómero reforzado con plata**
- PT** **Material de obturação de vidro ionómero reforçado com prata**
- NL** **Met zilver versterkt glasionomeer vulmateriaal**
- EL** **Υαλοϊονομερές κονία για εμφράξεις ενισχυμένη με άργυρο**
- SV** **Silverförstärkt glasionomer fyllningsmaterial**
- FI** **Hopeavahvisteinen lasi-ionomeeritayteenale**
- DA** **Sølvforstærket glasionomer-fyldningsmateriale**
- NO** **Sølvforsterket glassionomer fyllingsmateriale**

Instructions for Use
Gebruiksaanformatie
Mode d'emploi
Informazioni per l'uso
Información de empleo
Instruções de uso

**3M Deutschland GmbH**
Dental Products
Carl-Schurz-Str. 1
41453 Neuss – Germany

3M ESPE Customer Care/MSDS Information:
U.S.A. 1-800-634-2249 and Canada 1-888-363-3685.

3M, ESPE and Ketac are trademarks of 3M or 3M Deutschland GmbH.
Used under license in Canada. © 2014, 3M. All rights reserved.

- EN** **ENGLISH**

#### Description of Product

Ketac™ Silver is a silver-reinforced glass ionomer cement. On account of its chemical adhesion to the dental enamel and dentin, a gentle preparation and an especially tight seal at the margins is possible. The material can be applied without liting; it releases fluoride ions, and is radiopaque. The fillings are finishable to a faint lustre. The restorative material does not corrode and is physiologically as compatible as glass ionomer restorative material without silver.

⚠️ "For details on Ketac™ Conditioner and Ketac™ Glaze please refer to the respective Instructions for Use. The Instructions for Use must be kept for the duration of product use.

#### Indications

- Linings under amalgam fillings
- Core build-ups
- Uni-planar fillings and cervical fillings, if aesthetics is not the prime consideration
- Primary tooth fillings

#### Precautionary measures

3M MSDSs can be obtained from www.mmm.com or contact your local subsidiary.

#### Preparation

- Remove carious tooth substance only; undercuts are not required.
- Do not prepare thin bevelling margins. A wall thickness of 0.5 mm has to be secured for the material.

#### Pulp Protection

- In order to prevent pulpal irritation cover areas in close proximity to the pulp with local application of a calcium hydroxide material.

#### Conditioning

- For optimal chemical bonding to the teeth the smear layer produced in the preparatory steps must be carefully removed by applying Ketac Conditioner to the prepared surfaces and allowing the substance to react for 10 sec.
- Then rinse with copious amounts of water.
- Blow the cavity dry in only 2–3 short intervals with water and oil-free air or dab dry with cotton pellets. Do not over-dry! **The cavity should be just dry that the surface has a matt shiny appearance.** Excessive drying can result in post-operative sensitivity after filling.
- Avoid further contamination.

#### Dosing

The recommended mixing ratio by weight is 3.8 parts powder (approx. 1 spoonful) to 1 part liquid (1 drop). Please keep in mind that the enclosed dosing spoon is a guideline for dosing only. Add as much powder as possible to produce a ductile, packable consistency.

- Loosen the powder by shaking the bottle.
- Dip the enclosed spoon into the powder and use the plastic scraper at the bottle opening to level it off. Avoid compressing the powder.
- Dose a sufficient amount of powder on a mixing pad.
- Hold the bottle of liquid in a vertical position and dose the appropriate number of drops for the amount of powder next to the powder.
  - There must not be any dried residue of the liquid on the nozzle, and the drops must not contain any air bubbles!
- Carefully re-close both bottles after use and ensure that no powder residue remains on the plastic scraper of the powder bottle.

- Mixing
Ketac Silver at an ambient temperature of 20 °C to 25 °C/68 °F to 77 °F.
- Use a cement spatula made of metal or plastic on a mixing pad or a glass slab for mixing. Repeat the paste several times until a homogeneous consistency is obtained.

#### Application

- During the application and the setting phase, the working area should be protected from water and saliva. Avoid desiccation of the fresh filling, e. g. by a strong air stream.
- When using the material as a build-up under amalgam, the cavity has to be filled in such a way that the dentin is completely covered.
- Ketac Silver adheres to metal instruments, and for this reason should be washed off with cold water before setting.

#### Times

The following times apply for an ambient temperature of 23 °C/73 °F and a relative humidity of 50%:

Mixing	min:sec
	0:45
Working (including mixing)	4:00
Setting time (from the beginning of mix)	7:00
Setting time in the mouth	3:00

Higher temperatures shorten, lower temperature extend the working time (e.g. mixing on a cooled glass slab). A higher amount of powder also shortens the working time. An over-extended working time will cause the loss of adhesion to the dental enamel and the dentin.

#### Filling Protection

Glass ionomers set in two phases and are susceptible to moisture during the first setting phase. Subsequently dehydration may occur; which can be prevented by applying Ketac Glaze immediately after removal of strips or matrix.

- Apply Ketac Glaze to all exposed areas of the glass ionomer cement and light-cure for 10 sec.
- The inhibitory layer that remains after the Ketac Glaze has cured can prevent the curing of additional-type silicone impression materials. Either remove this layer or use a polyether impression material.

#### Finishing

- Remove excess with a sharp blade.
- Finishing should be carried out using rotating instruments 7 minutes after start of mixing at the earliest.
- For fine finishing, use discs or rubber cups and water-cooling.
- Requested, apply Ketac Glaze after surface sealing and light-cure for 10 sec.

#### Storage and Shelf-Life

Protect the powder from moisture.
Store the product at 15–25 °C/59–77 °F.
Do not use after the expiry date.

#### Customer Information

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

#### Warranty

3M Deutschland GmbH warrants this product will be free from defects in material and manufacture. 3M Deutschland GmbH MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M Deutschland GmbH's sole obligation shall be repair or replacement of the 3M Deutschland GmbH product.

#### Limitation of Liability

Except where prohibited by law, 3M Deutschland GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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